

Roscoe Bath Safety Commodes

Three-In-One Commode

- Versatile 3-in-1 design allows for use as a raised toilet seat, bedside commode or toilet safety frame
- Durable steel-welded frame is attractive, strong and easy to clean
- Tool-free, easy-to-assemble, includes plastic snap-on seat with lid
- Easily raised or lowered from 15 ½" to 23 ½" in 1" increments
- Plastic arms provide added comfort
- Comes complete with commode bucket with handle, cover and splash shield
- Color: white frame, gray commode
- Packaged 4 per case
- Limited Lifetime Warranty
- HCPCS Code: E0163*



Commode Accessories	Item #	Qty
Universal bucket with handle and cover, gray	BTH-BUC	12
Splash guard, gray	BTH-SPL	6
Toilet seat with lid, elongated, gray	BTH-SLD	4
Toilet seat with lid, round, gray	BTH-RSLD	4

Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Adjustable Seat Height	Weight Limit
BTH-31C (elongated seat)	17"	22"	13 ½"	16 ½"	15 ½" – 23 ½"	300 lbs.
BTH-RD31 (round seat)	17"	22"	13 ½"	15"	15 ½" – 23 ½"	300 lbs.

Three-In-One Folding Commode

- Versatile 3-in-1 design allows for use as a raised toilet seat, bedside commode or toilet safety frame
- Durable steel-welded frame is attractive, strong and easy to clean
- Folds flat for convenient storage and travel
- Sturdy plastic seat snaps on for easy tool free assembly
- Easily raised or lowered from 15 ½" to 23 ½" in 1" increments
- Comes complete with commode bucket with handle, cover and splash shield
- Color: white frame, gray commode
- Packaged 1 per case
- Limited Lifetime Warranty
- HCPCS Code: E0163*



Commode Accessories	Item #	Qty
Small bucket for folding commode, gray	BTH-SMBUC	12
Small splash guard for folding commode, gray	BTH-SMSPL	6

Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Adjustable Seat Height	Weight Limit
BTH-FLCOM	17"	22"	13 ½"	15"	15 ½" – 23 ½"	300 lbs.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com

Roscoe Bath Safety Commodes

Heavy Duty Commode

- The versatile 3-in-1 design allows for use as a raised toilet seat, bedside commode or toilet safety frame
- Durable steel-welded frame is attractive, strong and easy to clean
- Tool-free, easy-to-assemble, includes plastic snap-on seat with lid
- Easily raised or lowered from 17 ½" to 24 ½" in 1" increments
- Plastic arms provide added comfort
- Comes complete with commode bucket with handle, cover and splash shield
- Color: gray
- Packaged 2 per case
- Limited Lifetime Warranty
- HCPCS Code: E0168*



Commode Accessories	Item #	Qty
Universal bucket with handle and cover, gray	BTH-BUC	12
Splash guard, gray	BTH-SPL	6
Toilet seat with lid, elongated, gray	BTH-SLD	4

Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Adjustable Seat Height	Weight Limit
BTH-HV31	26"	31"	13 ¾"	16 ½"	17 ½" – 24 ½"	600 lbs.

24-inch Extra Wide Heavy Duty Commode

- Extra wide platform seat is comfortable and maximizes stability
- Commode opening is contoured and indented in the front to accommodate the needs of larger users
- Welded steel construction and 24-inch snap-on platform designed for strength and stability
- Comes complete with commode pail, cover and splash shield
- Can function as a free standing commode, as a toilet safety frame, or as a raised toilet seat with arms
- Height readily adjusts from 16 ½" to 23 ½" in 1" increments
- Color: gray
- Packaged 2 per case
- Limited Lifetime Warranty
- HCPCS Code: E0168*



Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Adjustable Seat Height	Weight Limit
BTH-COM	24"	26"	23"	18"	16 ½" – 23 ½"	400 lbs.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com

Roscoe Bath Safety Commodes

24-inch Heavy Duty Drop-Arm Commode



- Features a reinforced steel frame and an extra wide snap-on platform seat designed for comfort and stability
- Equipped with easy-to-release arm mechanism allowing for safe lateral patient transfers to and from the commode
- Comes complete with commode pail, cover and splash shield
- Height readily adjusts from 20" to 24" in 1" increments
- Color: gray
- Packaged 2 per case
- Limited Lifetime Warranty
- HCPCS Code: E0168*

Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Adjustable Seat Height	Weight Limit
BTH-DAC	23"	26 ½"	29"	18 ½"	20" – 24"	600 lbs.

26-inch Bariatric Drop-Arm Commode



- Features a reinforced steel frame and an extra wide padded platform seat designed for comfort and stability.
- Equipped with easy-to-release arm mechanism allowing for safe lateral patient transfers to and from the commode
- Plastic armrests provide added comfort and support
- Comes complete with commode bucket and cover
- Color: white
- Packaged 1 per case
- Limited Lifetime Warranty
- HCPCS Code: E0168*

Item Number	Width Between Arms	Outside Width	Seat Width	Seat Depth	Seat Height	Weight Limit
BTH-26BDAC	26"	29 ¼"	26"	21"	20 ½"	800 lbs.

Commode Accessories	Item #	Qty
Universal bucket with handle and cover, gray	BTH-BUC	12
Splash guard, gray	BTH-SPL	6
Toilet seat with lid, elongated, gray	BTH-SLD	4
Toilet seat with lid, round, gray	BTH-RSLD	4
Small bucket for folding commode, gray	BTH-SMBUC	12
Small splash guard for folding commode, gray	BTH-SMSPL	6

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com